1 Deutsch

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5 English

Description of this Thermometer

- Standby button ON/OFF / 2. LCD / 3. Measuring probe / 4. Protective cover

Instruction Manual

Read and save these instructions

Before using your device, please read this instructions and warnings completely. Remove the transparent protective layer from the thermometer panel before the first use.

Features

- Stimulation of nipple
The thermometer probe of the nipple adopts simulation design, accord with human body engineering.
- Safe and non-toxic
Use of medical and food grade silicone material, safe non-toxic. The device comes up to IP24, easy to clean.
- Easy to operate with accurate reading
Control by micro - computer with Buzz.
- Friendly environmental protection
No glass, no mercury, no pollution to environment and skin.
- Battery: No hazards for babies
For safety reasons the device cannot be opened, thereby preventing hazards of children swallowing the battery.

Safety Precautions

Marks or Symbol

	MUST OBSERVE Means Obligatory with detailed items expressed in words or figures within or beside the mark. Left one means General Compulsory.
	PROHIBITION Means Forbidden with detailed items expressed in words or figures within or beside the mark. Left one means General Forbithen.
	Means a possibility of personal injury in case of improper use.
	Means a possibility of personal injury or property damage in case of improper use.

Warning

Before the first use and after every successive use, clean and sanitize the thermometer's using medical alcohol.
• Use of thermometer without prior disinfection may communicate pathogenic bacteria.

This product is only used to measure the baby and children's oral temperature use, please do not use for other than the intended use. The use of expectations may cause injury or misleading information to you or patients.

For safely reasons, check the silicone for damage before and after each use.
• A damaged thermometer must not be used!

The thermometer is not a toy, store it in a safe place out of reach of children.
• Children may damage the thermometer.

Do not use the unexpected method of cleaning and disinfecting by cleaning and disinfection of thermometer.
• Unexpected cleaning and disinfection way may damage the equipment or cause fire.

The measuring result can only used for health monitoring and reference. Please don't conduct diagnose and treatment based on the measuring result, if you find the body temperature is abnormal, please contact your doctor.

Don't use the sharp tool, scissors etc. to scissor, cut, shear, scratch the measuring probe.

• When you use again, the measuring probe parts may be peeled off and swallowed by the bite and pull.

This product should not be placed in the near heat source or sources of ignition.
• High temperature may cause the product to the deformation of the casing, and even the fire.

6 English

Notice

Don't take temperature after drinking hot water.
• The hot water will enhance the oral temperature temporarily, it will cause the overhigh body temperature by immediately use.

Don't take the temperature immediately when you take the medicine which may enhance the body temperature.
The body temperature may be enhanced under the drug effect.

Don't talking and eating when you are measuring the oral temperature. And during the measure, use nose to breathe to prevent the measuring probe be exposed into air.

When open the mouth, the air will flow into and the probe is exposed which cause the inaccurate measuring result.

Don't take temperature after sports, bath or swimming.
The furious sports will enhance the body temperature. Please waiting for 15 to 30 minutes in the room, then to use the thermometer.

Protect the thermometer from impact and dropping!

The minimum measurement time until the beep is heard must be maintained without exception.

The operating ambient temperature range for this thermometer is 10-40 °C (50 - 104 °F).

The thermometer must be stored according to the technical specifications.

Strong electromagnetic fields can effect the correct functions of the thermometer.

Recommendations

- Don't use this thermometer for other purposes.
- It is forbidden to leave the product exposed to any chemical solvent, direct sunshine or high temperature.
- Do not measure while talking on the phone.
- Do not measuring while near the high frequency equipment (such as high-frequency electrotome, microwave ovens).

Intended Use

Except for measurement of oral temperature, don't use this product for other purposes. It can be used for anybody, e.g. for new-born and children without teeth.

How to Use the Thermometer

Display information	
5	6
1. Temperature value / 6. Unit of temperature (°C) / 7. Low battery symbol	
Err	1. Hardware failure 2. The ambient temperature exceeds the scope of design
Hi	The target value of measuring value is higher than 42.9 °C (109.9 °F)
Lo	The target value of measuring value is lower than 32.0 °C (89.6 °F)

Turning on the Thermometer

To turn on the thermometer, press the Standby/ Scan button (1), a short beep signals (thermometer ON). A display shelf-test is performed, all segments should be displayed. The test measurement reading will be shown on the display (2) automatically for 1 second. Then at an ambient temperature of less than 32.0 °C (89.6 °F), an „Lo“ and “°C“ appear at the display (2). The thermometer is now ready for use.

Function Test

Correct functioning of the thermometer is tested automatically each time it is turned on. If a malfunction is detected (thermistor is open-circuit or short-circuit), should be showing the "Err" in the LCD, and a measurement becomes impossible. In this case, the thermometer must be replaced.

High body temperature hint function

If the measurement value ≥ 37.8 °C (100.0 °F) a rapid beeping will be sound for 10 seconds to hint.

Using the thermometer

To ensure the accuracy of the measurement, turn on the thermometer at first. Then put the probe inside the mouth to keep close status.

When taking a measurement, the current temperature is